

AP4 Retrieve Release Notes

AP4 Retrieve 1.1

18 September 2014

Contents

ABOUT AP4 RETRIEVE	2
INSTALLATION.....	2
<i>System Requirements</i>	<i>2</i>
<i>To Install AP4 Retrieve</i>	<i>2</i>
REVISION HISTORY	2
<i>AP4 Retrieve 1.1.....</i>	<i>2</i>
<i>AP4 Retrieve 1.0.....</i>	<i>2</i>

About AP4 Retrieve

AP4 Retrieve is the software for downloading data from the AP4 Porometer.

Installation

System Requirements

To install and use AP4 Retrieve you will need:

- A PC running Windows XP, Vista, 7, 8 or a later.
- One free RS232 serial port, or USB-RS232 adapter.
- AP4 RS232 cable (supplied with the AP4).

To Install AP4 Retrieve

You need local administrator rights to install AP4 Retrieve 1.1 on computers with Windows 2000, XP and Vista. Please contact your system administrator for further information.

Select the **Install AP4 Retrieve 1.1** link in the Contents page of the Delta-T Software and Manuals CD and follow the on-screen instructions.

Alternatively, download AP4 Retrieve 1.1 from www.delta-t.co.uk and click **Run** when the download completes.

The installation program creates an **AP4 Retrieve** shortcut in the Start menu and on the desktop.

Revision History

AP4 Retrieve 1.1

AP4 Retrieve 1.1 includes the following changes:

1. COM port list now shows all available ports on the computer instead of just COM1 to COM8.
2. Now supports COM ports above COM9.
3. Data file selection now uses standard Windows save file dialog.
4. Serial port settings changed to match AP4 firmware 2.5.

AP4 Retrieve 1.0

AP4 Retrieve 1.0 is the first release of AP4 Retrieve.